

# Dr. Sonic 210 - Med 210



IDEALI PER FORNIRE ARIA COMPRESSA A 3 RIUNITI.  
 IDEAL TO DELIVER COMPRESSED AIR TO 3 DENTAL CHAIRS.

- Serbatoio da 50 litri, con protezione interna antiossidazione.
- Ventola ad alta portata convogliata per un facile raffreddamento.
- Interruttore on-off.
- Cavo con spina Plug-and-Play.
- Elettrovalvola "secure starting".
- Filtro uscita aria.
- Antivibranti serbatoio a ventosa per una perfetta installazione.
- Doppio manico per una facile movimentazione.
- Protezione termica del motore.
- Versione silenziosa "Dr Sonic" con carenatura insonorizzante ed essiccatore.
- Versione "Med" completabile con cabina fonoassorbente ed essiccatore.

- 50-lt tank, with inner anti-oxidant coating.
- High-capacity fan and forced air system for improved cooling.
- On-off switch.
- Cable with Plug-and-Play socket.
- "Secure starting" solenoid valve.
- Air outlet filter.
- Tank vibration pads, with suction cups for a safer installation.
- Double handle for easier handling.
- Motor thermal cut-out.
- "Dr Sonic" silenced version with soundproofing casing and drier.
- "Med" version completable with dryer and sound-proofed casing cover.



Dr. Sonic 210-50V-ES-2



Med 210-50V-2

ℓ	Codice Code	Prodotto Product	⚡ Vol/Ph	⚙️ W	⚙️ HP	🏠 Model	🏠 l/min. m³/h c.f.m.	🏠 bar p.s.i.	🏠 R.P.M.	🏠 n°	🏠 n°	🏠 dB(A)	🏠 L x D x H (cm)	🏠 kg lbs
50	871QS4A722	Dr. Sonic 210-50V-ES-2	400/3	1500	2	VKM210	180 10,8 6.3	8 116	1420	2-V	1	62	70 x 50 x 95	70 154.3
50	871QS4A604	Dr. Sonic 210-50V-ES-2M	230/1	1500	2	VKM210	180 10,8 6.3	8 116	1420	2-V	1	62	70 x 50 x 95	70 154.3
50	871QS1A722	Med 210-50V-2	400/3	1500	2	VKM210	180 10,8 6.3	8 116	1420	2-V	1	72	60 x 50 x 87	50,5 111.3
50	871QS1A604	Med 210-50V-2M	230/1	1500	2	VKM210	180 10,8 6.3	8 116	1420	2-V	1	72	60 x 50 x 87	50,5 111.3

Versioni con diverse tensioni e/o frequenze disponibili a richiesta.  
 Versions featuring different voltage and/or frequency ratings are available on request.