

# Dr. Sonic 420 - 640



IDEALI PER FORNIRE ARIA COMPRESSA DA 6 A 8 RIUNITI.  
 IDEAL TO DELIVER COMPRESSED AIR FROM 6 UP TO 8 DENTAL CHAIRS.

- Disponibili solo in versione silenziosa "Dr. Sonic", con carenatura insonorizzante ed essiccatore.
- Serbatoio da 90 litri, con protezione interna antiossidazione.
- Ventola ad alta portata convogliata per un facile raffreddamento.
- Elettrovalvola "secure starting".
- Filtro uscita aria.
- Antivibranti serbatoio a ventosa per una perfetta installazione.
- Doppio manico per una facile movimentazione.
- Protezione termica del motore.
- Apparecchiatura elettronica di comando a risparmio energetico per utilizzare al meglio i compressori.

- Only available in "Dr. Sonic" silenced version, with soundproofing casing and drier.
- 90-lt tank, with inner anti-oxidant coating.
- High-capacity fan and forced air system for improved cooling.
- "Secure starting" solenoid valve.
- Air outlet filter.
- Tank vibration pads, with suction cups for a safer installation.
- Double handle for easier handling.
- Motor thermal cut-out.
- Energy-saving electronic control equipment to use compressors in the best conditions.



Dr. Sonic 640-90F-ES-6

	Codice Code	Prodotto Product	Vol/Ph	W	HP	Model	l/min.	m³/h	c.f.m.	bar	p.s.i.	R.P.M.	n°	n°	dB(A)	L x D x H (cm)	kg	Lbs
90	873QS4A722	Dr. Sonic 420-90F-ES-4	400/3	3000	4	2-VKM210	360	21,6	12.7	8	116	1420	4-V	1	64	120 x 78 x 82	114	251
90	873QU4A722	Dr. Sonic 640-90F-ES-6	400/3	4500	6	2-VKM320	470	28,2	16.6	8	116	1420	4-V	1	64	120 x 78 x 82	118	260

Versioni con diverse tensioni e/o frequenze disponibili a richiesta.  
 Versions featuring different voltage and/or frequency ratings are available on request.