

Codice Code	Prodotto Product
KIT ACCESSORI ACCESSORIES KIT	
448458200V	Kit essiccatore MK 110-160 / MK 110-160 drier kit
448458300V	Kit essiccatore VKM 210-320 / VKM 210-320 drier kit
499QS0315V	Kit super-silenziamento SERIE "MED" / "MED" line super-silencing kit
048458010	Confezione sali essiccatore / Drier salt - 1 pack.

Kit essiccatore composto da: staffe, scambiatore di calore completo di elettroventola, separatore di condensa con scarico automatico, essiccatore ad adsorbimento, collegamenti elettrici e pneumatici. **Kit super-silenziamento** composto da: cabina componibile in materiale fonoassorbente, elettroventola di raffreddamento comandata da termostato.

Dryer kit made up of: brackets, heat exchanger fan, condensate separator with automatic drain, adsorption dryer, accessories for electrical and pneumatic connections. **Super-silencer kit** made up of: plastic shells with soundproofing lining, cooling fan controlled by thermostat.



499QS0315V



448458200V

Codice Code	Prodotto Product	L/min. m ³ /h c.f.m.			MAX bar p.s.i.		L x D x H (cm)	kg Lbs	
ESSICCATORE A REFRIGERAZIONE REFRIGERATION DRIER									
548500000	RD 6	600	36	21	14	203	37 x 52 x 48	25	55



548500000

Codice Code	Prodotto Product	L/min. m ³ /h c.f.m.			MAX bar p.s.i.		Ø x h (cm)	BSP
FILTRI FILTERS								
548103000	QFI 05	500	30	18	16	232	9 x 22	1/2"
548118000	PFI 05	500	30	18	16	232	9 x 22	1/2"
548148000	CFI 05	500	30	18	16	232	9 x 22	1/2"



548103000

Codice Code	Prodotto Product	MAX bar p.s.i.		BSP
SCARICATORI DI CONDENZA CONDENSATE DRAIN				
548302000	T1	16	232	1/8"
548306000	AUTO-DRAIN 950	16	232	3/8"



548302000



548306000

DISTRIBUTORE / DISTRIBUTOR

FINI è associata:
FINI is a member of:



Unione Nazionale Industrie Dentarie Italiane

Valori di rumorosità ± 3 dB(A) secondo norme PNEUROP/CAGI PN8NTC 2.3 (rilevati a 1 m). I valori di aria resa sono rilevati secondo ISO 1217.

Noise level - 3 dB(A) as PNEUROP/CAGI PN8NTC 2.3 (1 m). Free air delivery as for ISO 1217.

* MEDICAIR: 3 ANNI GARANZIA sul gruppo pompante, 12 MESI sugli altri componenti. AIRLAB: 12 mesi.
* MEDICAIR: 3 YEARS WARRANTY on the compressor pump unit, 12 MONTHS on the other components. AIRLAB: 12 months.



FINI S.p.A. Via Toscana, 21 - 40069 Zola Predosa - Bologna - Italy - Tel. +39 051 6168111
Fax: +39 051 752408 - +39 051 7459588 (Uff. Comm.le Italia) - +39 051 7417272 (Export Dept.)
www.finicompressors.com - info@finicompressors.com

www.finimedicaire.com - info@finimedicaire.com

